




## miRNA Agomirs & Antagomirs

**Agomirs & antagomirs** are small synthetic RNAs with modifications including 5' end-(phosphorothioate)<sub>2</sub>, 3' end-cholesterol, 3' end-(phosphorothioate)<sub>4</sub> and whole sequence 2'-methoxy group.

**RGBioTech** offers high-quality agomirs and antagomirs that can be used for either in vitro cell experiments or in vivo animal injections.

Item	Chemistry	Modification	Function
Agomir	Double-stranded RNA	Yes	Up-regulation of miRNA
Antagomir	Single-stranded RNA	Yes	Down-regulation of miRNA



-----For Research Use Only. Not For Use In Diagnostic Procedures.-----

