





miRNA Agomirs & Antagomirs

Agomirs & antagomirs are small synthetic RNAs with modifications including 5' end-(phosphorothioate)₂, 3' end-cholesterol, 3' end-(phosphorothioate)₄ and whole sequence 2'-methoxy group.

RGBioTech offers high-quality agomirs and antagomirs that can be used for either in vitro cell experiments or in vivo animal injections.

Item	Chemistry	Modification	Function
Agomir	Double-stranded RNA	Yes	Up-regulation of miRNA
Antagomir	Single-stranded RNA	Yes	Down-regulation of miRNA



-----For Research Use Only. Not For Use In Diagnostic Procedures.-----

